

HeelZup™ CVI

Reliable Above-The-Heart Leg Elevation

For conditions such as edema and chronic venous insufficiency, clinicians often recommend elevating the lower legs to a height that is above the heart when lying in a supine position. Until now, elevating the legs to this height was a challenging task often addressed with improvised approaches, using materials such as pillows, sofa cushions or blankets. These materials lack the necessary therapeutic qualities for adequate pressure redistribution and an acceptable level of comfort. Improvised approaches often result in poor compliance to leg elevation protocols due to their inconvenience, instability or discomfort. The HeelZup® CVI Leg Elevating Cushion is designed to safely and reliably elevate the lower legs above the level of the heart while offering several key benefits to the patient:

- *Elevates* the lower legs above the level of the heart, (ten inches at the distal end).
- *Eliminates* pressure on the Popliteal Fossa.
- *Eliminates* Pressure on the heels.
- *Envelopes* the thigh and calf in a soft bed of Visco Elastic Foam for optimal comfort and pressure redistribution.
- Available in one and two leg width versions and various sizes.
- No-Slip fabric on the bottom of the cushion prevents migration away from the patient.



Intensive Therapeutics, Inc.

Patent pending


Phone: 866.736.8484
Fax: 941.359.6939
www.heelzup.com




P.O. Box 488
Sarasota, FL
34230



The PopliZone™ contour design eliminates pressure on the Popliteal Fossa, greatly improving comfort and therapeutic efficacy. The heels are suspended and free of pressure, as seen below.



The HeelZup® CVI Cushion is covered with a 4-way stretch cover that is moisture proof and comfortable against the skin. A no-slip fabric is used on the bottom panel to prevent the cushion from migrating away from the patient.



Beneath the 4-way stretch cover lies a soft, comfortable bed of Visco-Elastic Foam, which enhances comfort and pressure redistribution.

Ordering Information:

Order Number Dual Leg Width	Order Number Single Leg Width	Description	Patient Height
HZCVI-SM	HZCVI-SM-I	HeelZup® CVI Cushion, Reusable Cover, Poplize™ Contour, Small	5'3" - Under
HZCVI-M	HZCVI-M-I	HeelZup® CVI Cushion, Reusable Cover, Poplize™ Contour, Medium	5'3" - 5'7"
HZCVI-LG	HZCVI-LGI	HeelZup® CVI Cushion, Reusable Cover, Poplize™ Contour, Large	5'8" - 6'
HZCVI-XLG	HZCVI-XLG-I	HeelZup® CVI Cushion, Reusable Cover, Poplize™ Contour, Extra Large	6'1" - Taller

Warranty:

This warranty applies only to defects in structure, workmanship and materials and only to damages arising from normal hospital usage on a non-prorated basis for a period of three (3) years for reusable models from the date of purchase. Damages arising from abnormal use such as those caused by needle punctures, burns, chemicals, negligent use or improper care or cleaning are excluded from coverage. Manufacturer will pay shipping charges incurred in connection with this warranty during the first year, thereafter, those charges are the facilities' or dealer's sole responsibility.