


Bari-Gel™

Bariatric Pressure Reducing Wheelchair Cushion

- Unique layering system of gel and visco-elastic foam provides pressure reduction and enhanced comfort for chair-bound patients weighing up to 850 pounds.
- The pressure reducing properties of the Bari-Gel™ Cushion will prevent and treat up to Stage II Pressure Ulcers.
- Includes four ties, which anchor the Bari-Gel™ Cushion to the wheelchair or other seating surface.
- The bottom is covered with a non-skid material, which, when combined with the ties will prevent the Bari-Gel™ from shifting and sliding upon ingress and egress.
- Includes a “domed” gel chamber, which spans the entire front portion of the cushion. This chamber is designed to prevent patients from sliding forward while sitting in the cushion.
- Includes a six (6) month limited, non-prorated warranty.



Intensive Therapeutics, Inc.

Phone: 866.736.8484
Fax: 410.203.9013
Web: www.heelzup.com



P.O. Box 488
Sarasota, FL
34230

Bari-Gel™

Bariatric Pressure Reducing Wheelchair Cushion

Sizes:

BG2618 26 x 18

BG2818 28 x 18

BG2620 26 x 20

BG2820 28 x 20

BG3020 30 x 20

BG2622 26 x 22

BG2822 28 x 22

BG3022 30 x 22

Warranty:

This warranty applies only to defects in structure, workmanship and materials and only to damages arising from normal hospital usage on a non-prorated basis for a period of six (6) months from the date of purchase. Damages arising from abnormal use such as those caused by needle punctures, burns, chemicals, negligent use or improper care or cleaning are excluded from coverage. Manufacturer will pay shipping charges incurred in connection with this warranty during the first three months, thereafter, those charges are the facilities' or dealer's sole responsibility.

Intensive Therapeutics, Inc.

Phone: 866.736.8484
Fax: 410.203.9013
Web: www.heelzup.com



P.O. Box 488
Sarasota, FL
34230