

# HeelZup®

## Therapeutic Heel Elevating Cushion with PopliZone™ Thigh Support

One of the leading challenges for the wound care community is heel pressure ulceration. Heels are naturally susceptible to pressure ulcers and even the most advanced support surfaces prove to be ineffective at prevention and treatment. Existing modalities of heel pressure ulcer prevention and treatment have many drawbacks. Heel “boots” require frequent removal to inspect the skin integrity of the lower leg. Pillows, which slide out of place and compress very easily are also a high maintenance and unreliable solution to the problem. The HeelZup™ Therapeutic Heel Elevating Cushion with PopliZone™ Thigh Support provides the most practical solution to this complex problem:

- Elevates the heels while supporting the leg from the thigh to the ankle.
- Eliminates pressure on the Popliteal Fossa.
- Raised side bolsters prevent legs from falling off of the sides of the cushion.
- Stays neatly in place and does not require constant maintenance and adjustment.
- Coated nylon cover reduces friction and shear, while the No-Slip bottom material helps the cushion to stay neatly in place. No need for constant maintenance and adjustment.
- The PopliZone™ Thigh Support can be retrofitted to previously purchased HeelZup® cushions that are not equipped with the PopliZone™ Thigh Support.
- Hinged seam between the PopliZone™ Thigh Support and HeelZup® Cushion allows the entire system to articulate.



**Intensive Therapeutics, Inc.**

Phone: 866.736.8484  
Fax: 941.359.6939  
www.heelzup.com



U.S. Patent # 6,634,045  
Other Patent Pending

**P.O. Box 488  
Sarasota, FL  
34230**



The HeelZup® Cushion with PopliZone™ Thigh Support is covered with a moisture proof, friction and shear reducing cover. An available retrofit kit enables owners of first generation HeelZup® cushions to upgrade to the PopliZone™ Thigh Support. The PopliZone™ Retrofit Kit reduces cost and eliminates unnecessary waste.



The HeelZup® Cushion and the PopliZone™ Thigh Support are joined together by a reinforced hinge style seam. This enables the cushions to bend easily when the bed articulates or when the patient changes position. When the cushion is not in use it can be folded to consume less space in storage.

## Ordering Information:

Order Number	Description	Patient Height	Patient Weight
TSH-13RPS	HeelZup™ Reusable PetiteSoft™, 13" Depth with PopliZone™ Thigh Support	Up to 5'5"	150 Pounds or Under
TSH-14R <b>STANDARD SIZE</b>	HeelZup™ Reusable, 14" Depth with PopliZone™ Thigh Support	5'6" - 6' Tall	Up to 300 Pounds
TSH-15R	HeelZup™ Reusable, 15" Depth with PopliZone™ Thigh Support	Over 6' Tall	Up to 300 Pounds
TSH-14RB	HeelZup™ Reusable Bariatric, 14" Depth with PopliZone™ Thigh Support	5'6" - 6' Tall	Over 300 Pounds
TSHxxR-KIT	HeelZup® Reusable Cover with PopliZone™ Thigh Support Kit	Specify Size	Specify Size

## Warranty:

This warranty applies only to defects in structure, workmanship and materials and only to damages arising from normal hospital usage on a non-prorated basis for a period of six (6) months for disposable, covered models and three (3) years for reusable models from the date of purchase. Damages arising from abnormal use such as those caused by needle punctures, burns, chemicals, negligent use or improper care or cleaning are excluded from coverage. Manufacturer will pay shipping charges incurred in connection with this warranty during the first year, thereafter, those charges are the facilities' or dealer's sole responsibility.