


No-Slip Wedge™

The No-Slip Wedge™ completes the ITI Patient Positioner Line by adding a highly reliable lateral turning solution. The No-Slip Wedge™ is made of a 1.85 pound density foam that is soft and comfortable, yet supportive enough to hold the torso in a turned position. A unique contour-cut pattern on the patient contact surface permits air circulation as it conforms to the shape of the patient. The No-Slip Wedge™ is an ideal way to enhance your pressure management and turning protocols:

- High Resiliency foam offers ideal pressure redistribution characteristics.
- Reliable way to maintain 30 degrees of lateral rotation.
- Contour cuts permit air circulation between the No-Slip Wedge™ and the user.
- Available in both single patient use and reusable models.
- Manufactured of high-density, anti-microbial, hypo-allergenic and latex-free foam.
- Reusable version is covered with a skin friendly, moisture proof friction and shear reducing cover with a panel of "No-Slip" fabric on the bottom.
- Infection Control friendly cover design features and Dartex cover material option.
- Available in both 20" Standard Length and 26" Extra Long Length at no additional charge.
- No-Slip fabric on the bottom of the reusable version prevents the wedge from migrating away from the patient.



Reusable No-Slip Wedge

U.S. Patent # 7,240,384



Single Patient Use No-Slip Wedge

Phone: 866.736.8484
Fax: 941.359.6939
www.heelzup.com

Intensive Therapeutics, Inc.


P.O. Box 488
Sarasota, FL
34230

Specifications:



No-Slip Wedge Disposable:

Manufactured of high-density, anti-microbial, hypoallergenic foam. The dimensions are 20" W x 11" Deep x 7" Tall. The No-Slip Wedge™ meets or exceeds the specifications in California Upholstery Provisions as stated in Technical Bulletin #117.



No-Slip Wedge Extra-Long Disposable:

Manufactured of high-density, anti-microbial, hypoallergenic foam. The dimensions are 26" W x 11" Deep x 7" Tall. The No-Slip Wedge™ meets or exceeds the specifications in California Upholstery Provisions as stated in Technical Bulletin #117.



No-Slip Wedge Reusable:

Manufactured of high-density, anti-microbial, hypoallergenic foam. The dimensions are 20" W x 11" Deep x 7" Tall. Upholstered with a moisture proof, reusable, coated nylon fabric cover, which is friction and shear reducing. No-Slip fabric is used on the bottom to prevent migration away from the patient. The No-Slip Wedge™ meets or exceeds the specifications in California Upholstery Provisions as stated in Technical Bulletin #117.



No-Slip Wedge Extra-Long Reusable:

Manufactured of high-density, anti-microbial, hypoallergenic foam. The dimensions are 26" W x 11" Deep x 7" Tall. Upholstered with a moisture proof, reusable, coated nylon fabric cover, which is friction and shear reducing. No-Slip fabric is used on the bottom to prevent migration away from the patient. The No-Slip Wedge™ meets or exceeds the specifications in California Upholstery Provisions as stated in Technical Bulletin #117.

Infection Control Cover Features:



Waterfall Flap covers the top half of the zipper enclosure.



Raised border behind the bottom zipper prevents fluid infiltration.



Heat sealed seams available on optional Dartex® cover.

Ordering Information:

NSW-D	No-Slip Wedge, Disposable, Bare Foam
NSW-R	No-Slip Wedge, Reusable, Nylon Cover, No-Slip Bottom
NSW-DXL	No-Slip Wedge, Disposable Extra Long, (26" length), Bare Foam
NSW-RXL	No-Slip Wedge, Reusable Extra Long, (26" length), Nylon Cover, No-Slip Bottom



Warranty:

This warranty applies only to defects in structure, workmanship and materials and only to damages arising from normal hospital usage on a non-prorated basis for a period of six (6) months for disposable, covered models and three (3) years for reusable models from the date of purchase. Damages arising from abnormal use such as those caused by needle punctures, burns, chemicals, negligent use or improper care or cleaning are excluded from coverage. Manufacturer will pay shipping charges incurred in connection with this warranty during the first year, thereafter, those charges are the facilities' or dealer's sole responsibility.