

# ImmersaFlex™ Series 03 / Series 04 Wheelchair Seating for Skin Protection

With Genuine TEMPUR® Technology



## New from Intensive Therapeutics!

ImmersaFlex™ therapeutic seating cushions are designed to provide superior performance at an attractive price. ImmersaFlex™ Series 03 and 04 wheelchair cushions are made with an advanced formula of genuine TEMPUR® material, which is ideal for therapeutic seating applications. Genuine TEMPUR® material, combined with our unique ImmersaFlex™ foundation, provides the user with a fluid and dynamic feel while maintaining excellent postural stability. The superior immersion characteristics of the TEMPUR® material used in ImmersaFlex™ wheelchair cushions assures excellent pressure redistribution and skin protection.

### Industry Leading Features:

- Multi-density foundation provides firm stability for the thighs while maintaining optimal comfort and pressure redistribution under the ischium.
- Genuine TEMPUR™ Material by Tempur-Pedic Medical top layer provides superior pressure redistribution and comfort while maintaining excellent postural stability.
- Recovery Knit fabric cover stretches in all directions and has a billowy fit to permit excellent immersion.
- No-Slip bottom fabric and quick release locking safety belt keep ImmersaFlex™ Cushions securely attached to the wheelchair.
- Integrity 30™ Flame Retardant Liner exceeds Medicare flammability requirements while maintaining a contour conforming fit.
- Tested by an independent laboratory for immersion, interface pressure and durability over a simulated period of five years.
- Meets or exceeds the criteria for HCPCS Code E2603, (widths up to 20") and E2604, (widths > 22").



*The ImmersaFlex Series 03 Cushion from the back, showing the quick release buckle and from the front, showing the comfort enhancing waterfall seam.*



*An exploded view of the Series 03 / Series 04 cushion shows the top layer of genuine TEMPUR® material, the softer ischial immersion zone and it's surrounding firm foundation area for the thighs.*

## Intensive Therapeutics, Inc.

P.O. Box 488  
Sarasota, FL 34230  
866.736.8484

[www.ImmersaFlex.com](http://www.ImmersaFlex.com)  
[info@ImmersaFlex.com](mailto:info@ImmersaFlex.com)

*All ImmersaFlex™ Cushions are equipped with two features for preventing movement of the cushion during use. A high friction no-slip bottom fabric and a durable, quick release locking safety belt, pictured here.*





ImmersaFlex™ Cushion covers are fitted with a large wrap-around zipper enclosure, which enables easy removal and insertion of the cushion before and after laundering.



Waterfall seam on the front of the cushion cover eliminates contact between the thighs and cover seams.

## Sizes and Ordering Information:

### ImmersaFlex™ Series 03 Size Chart

Cushion Name	Stock Number	Size, (W x D x Thickness)	Suggested HCPCS Code
ImmersaFlex™ Series 03	IFX-03-1616-3	16 x 16 x 3	E2603
ImmersaFlex™ Series 03	IFX-03-1618-3	16 x 18 x 3	E2603
ImmersaFlex™ Series 03	IFX-03-1816-3	18 x 16 x 3	E2603
ImmersaFlex™ Series 03	IFX-03-1818-3	18 x 18 x 3	E2603
ImmersaFlex™ Series 03	IFX-03-2016-3	20 x 16 x 3	E2603
ImmersaFlex™ Series 03	IFX-03-1616-4	16 x 16 x 4	E2603
ImmersaFlex™ Series 03	IFX-03-1618-4	16 x 18 x 4	E2603
ImmersaFlex™ Series 03	IFX-03-1816-4	18 x 16 x 4	E2603
ImmersaFlex™ Series 03	IFX-03-1818-4	18 x 18 x 4	E2603
ImmersaFlex™ Series 03	IFX-03-2016-4	20 x 16 x 4	E2603
ImmersaFlex™ Series 03	IFX-03-2018-4	20 x 18 x 4	E2603

### ImmersaFlex™ Series 04 Size Chart

Cushion Name	Stock Number	Size, (W x D x Thickness)	Suggested HCPCS Code
ImmersaFlex™ Series 04	IFX-04-2216-3	22 x 16 x 3	E2604
ImmersaFlex™ Series 04	IFX-04-2218-3	22 x 18 x 3	E2604
ImmersaFlex™ Series 04	IFX-04-2220-3	22 x 20 x 3	E2604
ImmersaFlex™ Series 04	IFX-04-2416-3	24 x 16 x 3	E2604
ImmersaFlex™ Series 04	IFX-04-2418-3	24 x 18 x 3	E2604
ImmersaFlex™ Series 04	IFX-04-2420-3	24 x 20 x 3	E2604
ImmersaFlex™ Series 04	IFX-04-2216-4	22 x 16 x 4	E2604
ImmersaFlex™ Series 04	IFX-04-2218-4	22 x 18 x 4	E2604
ImmersaFlex™ Series 04	IFX-04-2220-4	22 x 20 x 4	E2604
ImmersaFlex™ Series 04	IFX-04-2416-4	24 x 16 x 4	E2604
ImmersaFlex™ Series 04	IFX-04-2418-4	24 x 18 x 4	E2604
ImmersaFlex™ Series 04	IFX-04-2420-4	24 x 20 x 4	E2604

**Series 03** cushions are reimbursed under HCPCS Code **E2603** and are manufactured in widths of up to 20". **Series 04** cushions are reimbursed under HCPCS Code **E2604** and are manufactured in widths of 22" and wider.